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The screenshot displays the software interface for temperature control. At the top, there is a navigation bar with 'Home', 'Edit', 'Data', and 'Help' buttons, along with a 'connected' status indicator. The main area features a 'Temperature Waveform' graph showing Temperature (°C) on the y-axis (0 to 20) and Time on the x-axis (0 to 39). The current temperature is 22.8 °C. Below the graph, there are control panels for 'Set Temperature', 'Segment Mode', and 'Formulation Mode'. The 'Segment Mode' panel is active, showing 'Mode' set to 'time' (selected from a dropdown menu), 'Target temperature' set to 0, and 'Mode property' set to 0. The 'Formulation Mode' panel shows 'Formula' set to 'test'. There are 'Start' and 'Stop' buttons, as well as 'Re-rec' and 'Refresh' buttons.

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
The interface consists of a top navigation bar with 'Home', 'Edit', 'Data', and 'Help' buttons, and a status bar on the right. The main area features a 'Temperature Waveform' graph with 'Temperature(°C)' on the y-axis (27.25 to 30) and 'Time' on the x-axis (0 to 27). Below the graph are three mode selection buttons: 'Set Temperature' (checked), 'Segment Mode', and 'Formulation Mode'. The 'Set Temperature' section includes a 'Temperature' input field (30), 'Target temperature' (0), and 'Mode property' (0). The 'Formulation Mode' section includes a 'Formula' dropdown (test). Control buttons for 'Start' and 'Stop' are also present.

First Screenshot: The status bar shows 'fixed temperature start operation'. The graph displays a rising temperature curve. The 'Re-rec' button is highlighted with a red box.

Second Screenshot: A dialog box appears with the text 'Do you want to re-recording?' and 'Yes' and 'No' buttons. The 'Yes' button is highlighted with a red box.

Third Screenshot: The status bar shows 'controller closed'. The graph is empty, and the temperature on the right is 29.3 °C.

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